

More freedom.
More confidence.
With mylife.



With its mylife Diabetescare brand, Ypsomed offers a comprehensive portfolio of products and services for people with diabetes. This allows users easy, discreet and reliable self-treatment. With mylife Diabetescare, self-treatment becomes simpler and easier.

Ypsomed Ltd. // 1 Park Court // Riccall Road // Escrick // North Yorkshire // YO19 6ED // info@ypsomed.co.uk // www.mylife-diabetescare.co.uk // Customer Care: UK 0344 856 7820 or ROI 1800 911 522 (Calls to this number are free of charge from mobile phones and landlines, plus your phone company's access charge.)

Abbr: TIRp: Time in range in pregnancy // 1. Lee TTM, et al. Automated Insulin Delivery in Women with Pregnancy Complicated by Type 1 Diabetes. N Engl J Med. 2023 Oct 26;389(17):1566-1578. // Disclaimer: The product images are for illustrative purposes only. // mylife and YpsocPump are registered trademarks of Ypsomed AG. // CamAPS is a registered trademark of CamDiab Ltd. // The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott and used with permission. // Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. in the United States and/or other countries. For safety information on mentioned products, see Instructions for Use or www.mylife-diabetescare.com/safety



mylife CamAPS FX

For pregnancy and pre-pregnancy.

The only AID algorithm designed and licenced for use in pregnancy. Proven to improve the health of pregnant people with type 1 diabetes and their babies.

mylife Loop powered by CamAPS FX

The only AID algorithm licenced for use in pregnancy

The mylife CamAPS FX app is personalised and adaptive; it constantly learns and quickly adapts to ever changing insulin needs. It adjusts to diurnal and day-to-day variations in insulin requirements and compensates for over and under-dosing of the meal bolus. For increased privacy and convenience, you can deliver a meal bolus directly from your smartphone.



- Glucose levels are measured by the Dexcom G6 or Abbott FreeStyle Libre 3 continuous glucose monitoring systems.
- mylife CamAPS FX adjusts insulin delivery level every 8 to 12 minutes, based on the current and the predicted glucose levels.
- Insulin is delivered with the mylife YpsoPump.

Clinical results

AiDAPT randomised control trial¹

- Using mylife CamAPS FX during pregnancy improved glycaemic control:
 - 10.5% more time spent in the pregnancy specific range of 3.5 – 7.8 mmol/l (63-140mg/dL).¹
 - 0.3% reduced HbA1c.¹
 - No difference in time spent in hypoglycaemia.¹
 - The women in the closed loop group achieved the pregnancy HbA1c target of 6.0% (42 mmol/mol).¹
- These improvements were visible immediately after initiation of the mylife CamAPS FX and were sustained throughout the pregnancy.¹
- All reported more enjoyable pregnancy experiences as a result of using hybrid closed loop system.¹

“

“In terms of enjoyment factor, and how I feel about the baby, and how I feel about the coming labour, and how excited I am to meet them and the bonding process and stuff like that, I'd say I've had a lot more time for it this time round”¹

Participant 016

Find out more here

Scan the QR code to find out more



Website

Scan here to find out more about our mylife Loop system on our website.



Video

Scan here to watch a video on managing pregnancy with the mylife CamAPS FX app.



Clinical trials

Scan here to download the results of our clinical trials.